


ACCUGENCE® PLUS



Multi-Monitoring System

Model PM800



-IN-1 testing



ACCUGENCE® PLUS Multi-Monitoring System (Model No. PM 800) is an easy and reliable Point-of-Care meter which is available for Blood Glucose (GOD and GDH-FAD enzyme both), β -ketone, uric acid, hemoglobin testing from whole blood sample for the hospital primary care patients self-monitoring.

- + Multiple parameter measurement for Blood Glucose, β -ketone, uric acid, hemoglobin testing
- + User friendly design and reliable result with clinically proven accuracy
- + Glucose result range indicator and multi-function for better management

ACCUGENCE® PLUS

**Reliable result:**

Clinically proven accuracy result with plasma-equivalent calibration

**Fast testing time****User friendly:**

Less pain with tiny blood samples, allow redoes blood

**Advanced features:**

Before/after meal markers, hypo / hyper and ketone alarms, 5 daily testing reminders

**Intelligent identification:**

Intelligent recognize the test strips type, samples type or control solution

**Smooth strips ejector****Glucose result range indicator****Technical Specifications**

Product name	ACCUGENCE® PLUS Multi-Monitoring System			
Model No.	PM800			
Parameter	Blood Glucose (GOD & GDH), β -Ketone (KET), Uric Acid (UA), Hemoglobin(Hb)			
Measurement Range	GLU: 0.6 ~ 33.3 mmol/L (10 ~ 600mg/dl) KET: 0.0 ~ 8.0 mmol/L UA: 3.0 ~ 20.0 mg/dL (179 ~ 1190 μ mol/L) HB: 3.0 ~ 26.0 g/dL (1.9 ~ 16.1 mmol/L)			
Hematocrit Range	GLU and KET: 10% ~ 70 % UA: 25% ~ 60% HB: 9% ~ 77%			
Result Calibration	Plasma-equivalent			
Sample	GDH, KET, UA, HB: fresh capillary whole blood and venous Blood GOD: fresh capillary whole blood only			
Sample Size	GLU: 0.7 μ L KET: 0.9 μ L UA: 1.0 μ L HB: 1.2 μ L			
Test Time	GLU: 5s KET: 5s UA: 15s HB: 15s			
Memory	GLU (GOD + GDH): 300 tests KET: 100 tests UA: 100 tests HB: 200 tests			
Automatic Shutoff	2 minutes			
On/Off Source	Two CR 2032 3.0V coin cell batteries			
Battery Life	Around 2000 tests			
Operating Temperature	GLU & KET: 5 ~ 45 °C (41 ~ 113 °F) UA & HB: 10 ~ 40 °C (50 ~ 104 °F)			
Operating Humidity	10 ~ 90% (non-condensing)			

**e-LinkCare Meditech Co.,Ltd.**

Unit 46 Cariocca Business Park 2, Sawley Road Miles Platting, Manchester,
United Kingdom M40 8BB
E-mail: info@e-linkcare.com



e-LinkCare Meditech Co.,Ltd.
No. 30 Baita Tongjiang Road,
Taizhou, Zhejiang, China,
317317

www.e-linkcare.com